

ASSISTIVE LISTENING, ALERTING & TELEPHONE TECHNOLOGY IDEAS


ESTHER KELLY
HEARING LOSS RESOURCE SPECIALIST
esther.kelly@dactexas.org



**Funded by: Texas D.A.R.S.-
Office for Deaf and Hard of Hearing Services**



LISTENING TECHNOLOGY






T-Coils

USES OF T-COIL- when t-coil switch is activated:

- Better telephone listening when phones are t-coil compatible.
- Cell phones with induction loops
- Assistive technology such as FM and Infrared Listening Systems.
- T-Coils should be placed on every hearing aid so that this technology is available as current technology and future technology can be used in the work place, home and social environments.



► A dime with T-Coils

FM Systems

Comfort Contego

- ▶ Dual microphone
- ▶ More control of feature selections
- ▶ Omni and directional microphones



Office Phones ...New Ideas!

Try telephone adaptor from Radio Shack :

- **Comfort Contego** – plus Radio Shack adaptor
 - Part # 43-1237



Attached to Office Digital Phone Lines



Williams Sound Pocketalker Ultra



Serene DT-100

DirecTalk Hands-free Personal Listener



Clipboard “Carry Around” Hearing Loop



VIVID Acoustics V-A 1303

Loop Hear

By Geemarc



E- Scope

Electronic Stethoscope



E-Scope

Stethoscope Options-

Stethoscope Use with BTE Hearing Aids & Cochlear Implants



Ear Protection with Hearing Enhancement



**Hearing Protection
Earmuffs**



**Walker's Electronic
Power Muffs**



**Williams Sound
HED008**


ClearSounds® CS-CLA7v2



Professional Office Neckloop Systems

ClearSounds®

CONNECTIONS



- Multi-Line VoIP Office Phone
- Plantronics M22 Headset Base
- Microphone
- CLA7v2 Neckloop
- 30dB Adjustable Volume Control
- M22 Quick Connect Cable – Sold Separately.

\$249.95

BLUETOOTH TECHNOLOGY

**CLEAR
SOUNDS
Quattro
CLA7BT**



We can do what they do ...



Cellular



TV Listening



Music

... at a fraction of the cost.



ClearSounds®
Hear how good life can be

Call toll-free 800.965.9043 today!



ClearSounds
QUATTRO



RSO62 Earbuds



Bluetooth Dongle
for TV

ClearSounds QLink



Amplified Bluetooth Transmitter



Another Option for Office Phones



Univox 2A+ Chair /Desktop Loop




Serene VM-150



Amplified answering machine with slow play-back



Alerting Technology



Gentex Smoke Detector

Photoelectric Smoke Alarm with Visual Signaling Appliance



Nutone LA204WH

Wireless Doorbell and Receiver signals a ring at the doorbell



© Copyright 2010 Tellex, Inc.



Vibrating Alarm Clocks



Serene VA3

Ameriphone Big Time Alarm Clock



Sonic Ring Elite Telephone Signaler



Ringmaster Super Loud Cell Phone Ringer/Flasher Model RF110



Alerts you to telephone ringing and
text messages



Other Resources





Phone Communication Options

Web-CapTel Now Available

www.hamiltoncaptel.com

www.sprintcaptel.com





Answer the call
—with confidence.

<https://www.captioncall.com>





ClearCaptions

Apps for cell phone

<http://www.clearcaptions.com>





www.captionfish.com

OC - *Open Captioned* - The movie's text is superimposed over the screen images.

RW - *Rear Window[®] Captioned* - A transparent acrylic panel attached to your seat reflects the captions so that they appear superimposed on the screen.

USL - *USL Closed Captioned System (CCS)[®]* - The CCS is designed to enhance the hearing impaired cinema patron's movie-going experience. A single infrared emitter broadcasts closed caption text and two channels of audio into an auditorium which can be picked up by either a display that can be flexibly positioned in front of you or special eyewear.

CV - *CaptiView[®] Closed Captioning* - The CaptiView system consists of a small OLED display on a bendable support arm that fits into the theater seat cup holder. The easy-to-read screen is equipped with a high contrast display that comes with a privacy visor so it can be positioned directly in front the movie patron with minimal impact or distraction to neighboring patrons.

ST - *Subtitled* - Textual versions of the dialogue are displayed in English on the bottom of the screen.





CART: Communication Access Real-time Translation

- Communication service for persons with hearing loss provided by a specially trained realtime court reporter.
- Persons who are late-deafened, hard of hearing or deaf and have not learned sign language
- CART reporters write the spoken word.
- CART is used in classrooms, medical settings, churches, the work place and for meetings.
- The Americans with Disability Act includes CART in its definition as an “auxiliary aid”. This means one has the same rights to CART as they would to a sign language interpreter or Assistive Listening device.

TEXAS COURT REPORTERS ASSOCIATION

PO Box 40429

Austin, Texas 78704-0429

512 462-3650



Technology Catalogs

Teltex, Inc. 1 888-515-8120
1081 West Innovation Dr.
Kearney, MO 64060

Harris Communications 1 800-825-6758
15155 Technology Drive
Eden Prairie, NM 55344

Hearing Systems Inc 1 877 855-8916
16103 W. Little York, Suite F
Houston, TX 77084

LS&S Products 1 800-468-4789
PO Box 673
Northbrook, IL 60065

Hitec 1 800 288-8303
8160 Madison Ave
Burr Ridge, IL 60527

Becky Morris 1 800 838-1649
Beyond Hearing Aids, Inc
463 Erlanger Rd., Suite 1
Erlanger, KY 41018



Questions?

Esther Kelly

Hearing Loss Resource Specialist

esther.kelly@dactexas.org

