

HEAR HERE

This symbol indicates that a **HEARING LOOP** is available.

Activate the telecoil (t-coil) in your hearing device to hear sound with less background noise — directly into your device.

If you do not have a t-coil enabled device, request a headset + receiver.

For questions about hearing loops and hearing loss contact or visit the Hearing Loss Association of America:

GITHLinfo@hearingloss.org
www.hearingloss.org



Hearing Loop
Switch hearing aid to T-coil

A LIFE-CHANGING SOLUTION FOR PEOPLE WITH HEARING LOSS

HEARING LOOPS

- **Simple** – hearing loops reduce background noise and send sound from a microphone directly to t-coil enabled devices
- **Discreet** – there is no need to ask for additional equipment
- **Effortless** – most hearing devices have an external toggle switch to enable the t-coil

Talk to your audiologist about t-coil equipped hearing devices. If you already have one of these devices, ask if it is t-coil enabled and learn how to activate the t-coil program.



Microphone



Hearing Aids
Cochlear Implants
Bone Conduction

Toggle to access the t-coil